

Hycell 2x qPCR Master Mix

特色

- ✓ **Compatible for various realtime cyclers**
更多適用機型
- ✓ **Highly Specific** 有效降低 primer dimer
- ✓ **Efficient** 良好的 PCR 效率
- ✓ **Broad Dynamic Range** 有效定量範圍更大
- ✓ **More Sensitive** 敏感度提高
- ✓ **SYBR green based、Probe based qPCR**
- ✓ 可以根據使用機型添加 ROX dye

5.04ml 只要 \$7,800

買1送1


適用於以下各種機型

品名	適用機型	
2X qPCR master mix with ROX	Applied Biosystems	ABI Prism 7000
		ABI Prism 7300
		ABI Prism 7700
		ABI Prism 7900
		StepOne Plus
		StepOne™
	Eppendorf	Mastercycler ep realplex4
2x qPCR master mix with low ROX	Applied Biosystems	ABI Prism 7500
	Stratagen	Mx3000
		Mx3005
		Mx4000


品名	適用機型	
2x qPCR master mix	BioRad	CFX96
	Roche	LightCycler 480
	MJ Research	Opticon
		Opticon 2
		Chromo 4
	Carbett Life Science	Roter-gene 600
		Roter-gene 3000
	Eppendorf	Mastercycler ep realplex2

靈敏度提高、定量準確度提高、再現性提高

PCR sensitive comparison



Quantitative qPCR for human β actin



Human β -actin mRNA was detected with serially diluted cDNA from HeLa cell total RNA.

另有 Tag Man probe 合成服務特價中，歡迎來電洽詢



翰新國際有限公司
Hycell International Co., Ltd.

TEL: 02-28239957
FAX: 02-28236645