

# GOAL Bio Equipments

## 實驗儀器

### 電動吸取器

Goalbio pipette (Pipet aid)



安全保護裝置  
protective hole



無段變速功能  
Continuous speed control wheel



可站立於實驗桌  
Self-inverted on the bench

Hydrophobic  
0.2 um filter

Helix Switch  
Closer-knit

Nose Cone

Battery Lidicator

Lithium Battery

Speed control wheel

Blow-out/gravity  
Dispense switch

大容量鋰電池，無記憶效應  
Lithium battery, Nomemory effect

DC-7u Socket

重量：167g，超輕量機種！  
Light-weighted: 167g

特惠價 \$7,800

再送 10 ml 刻度吸管

可抗紫外線照射  
UV-resistant housing

### 桌上型離心機

Cubee Mini Centrifuge



特惠價 \$4,000

再送1.5ml 離心管



專利三合一轉子  
Patented 3-in-1 rotor design

2~1.5ml tube x 8孔  
0.5ml tube x 8孔  
0.2ml tube x 16孔

橘  
Orange

粉紅  
Pink

粉藍  
Blue

▶ 使用4700rpm高轉速馬達，低噪音  
4,700 rpm speed with noise reduction design.

▶ 強化上蓋支撐點，耐用耐操  
Reinforced hinge for long-term durability.

▶ 底部加重，運作不位移  
Stabilized base to avoid movement during operating.

▶ 單手即可操作，簡單方便  
Open/Close design, capable for one-hand operation, safe and easy to operate

GOAL Bio

www.goalbio.com

HC 翰新國際有限公司

HyCell International Co., Ltd.

TEL: +886-2-2823-9957 FAX: +886-2-2823-6645