

Thermo
SCIENTIFIC

Haraeus™ Pico™ / Fresco™ 桌上型微量(冷凍)離心機



	Pico™ / Fresco™ 17	Pico™ / Fresco™ 21
Max RCF	17,000 xg	21,100 xg
Max Speed	13,300 rpm	14,800 rpm
Time Set Range	1 - 99 min; 1 min increments	
Temperature Range	-9 to +40°C; 1°C increments (Fresco only)	

極輕巧 · 極強大
體積小，效能高，節省寶貴實驗空間

生物安全防護
ClickSeal™ 上蓋有效杜絕汙染




人體工學設計
One-Click 開啟轉子及開關蓋
直觀極簡操作
自動轉子識別及不平衡檢測



24 x 1.5/2.0 mL



36 x 0.5 mL



18 x 2.0 & 0.5 mL



PCR 4 x 8



PCR 8 x 8



Hematocrit

ClickSeal™ Bio-Containment Lid

Screw-On Lid

▶ 六種不同轉子，日常實驗需求一機滿足！

HC 翰新國際有限公司



TEL: (02) 2877-1122
hycellbio@gmail.com



Haraeus™ Megafuge™

桌上型多功能離心機

Fits in.

- ▶ 機身纖巧，適合所有實驗室環境
- ▶ 相容性高，滿足您任何研究需求




容量大
極方便

Auto-Lock™ 專利
一鍵快速更換轉子



耐用·極輕-碳纖維
Fiberlite™ 轉子



高轉速



各式離心、採血管
96孔盤等皆可操作



專利ClickSeal™
生物防護蓋杜絕汙染

安全

Stands out.

- ▶ 獨特高容量轉子，節省重複離心之實驗時間
- ▶ SMARTSpin® 轉子自動辨識及加減速最佳化
- ▶ 直覺化操作及程式編輯，一鍵呼叫常用程式




Megafuge 8/8R

▶ 可離心96孔盤之最小機種



Megafuge 16/16R

▶ 1.6L容量，適合任何日常用途



Megafuge 40/40R

▶ 市面上最大容量(4L)之桌上型機型

