



陪您一起遨遊實驗世界

## 高品質二級抗體現正熱賣中！

Fluorescent	Color	EX wavelength (nm)	EM wavelength (nm)	Fluorescent groups
Dylight 405	Blue	400	421	Alexa 405 , AMCA
Dylight 488	Green	493	518	Alexa 488, FITC
Dylight 549	Yellow	550	568	Alexa 546, Alexa 555, TRITC, Cy3
Dylight 594	Orange	591	616	Alexa 594, Texas Red
Dylight 649	Red	652	670	Alexa 647, Cy5
Dylight 680	Near Red	682	715	Alexa 680, Cy5.5
Dylight 800	Infrared	770	794	IRDye 800

多種螢光及紅外光標的，滿足您所有實驗需求！

Cat. No.	Product Name	Cat. No.	Product Name	Cat. No.	Product Name
BS13271	Goat Anti-Rabbit IgG (H+L) Unconj.	BS50950	Goat Anti-Mouse IgG (H+L) FITC	BS50650	Goat Anti-Mouse IgG (H+L) Biotin
BS13278	Goat Anti-Rabbit IgG (H+L) HRP	BS30350	Rabbit Anti-Sheep IgG (H+L) HRP	BS10042	Rabbit Anti-Goat IgG(H+L) Dylight 800
BS50350	Goat Anti-Mouse IgG (H+L)-HRP	BS10008	Rabbit Anti-Goat IgG(H+L) Cy3	BS10039	Rabbit Anti-Goat IgG(H+L) Dylight 680
BS12471	Goat Anti-Mouse IgG (H+L) Unconj.	BS10032	Goat Anti-Rat IgG(H+L) Dylight 594	BS10031	Rabbit Anti-Goat IgG(H+L) Dylight 594
BS12478	Goat Anti-Mouse IgG (H+L) HRP	BS10025	Goat Anti-Rat IgG(H+L) Dylight 549	BS10024	Rabbit Anti-Goat IgG(H+L) Dylight 549
BS30503	Rabbit Anti-Goat IgG (H+L) HRP	BS10021	Goat Anti-Rat IgG(H+L) Dylight 488	BS10019	Rabbit Anti-Goat IgG(H+L) Dylight 488
BS10550	Goat Anti-Rabbit IgG (H+L) AP	BS10014	Goat Anti-Rat IgG(H+L) Dylight 405	BS10013	Rabbit Anti-Goat IgG(H+L) Dylight 405
BS50550	Goat Anti-Mouse IgG (H+L) AP	BS10036	Goat Anti-Rat IgG (H+L) Dylight 649	BS10035	Rabbit Anti-Goat IgG (H+L) Dylight 649
BS10001	Goat Anti-Rat IgG(H+L) TRITC	BS10041	Goat Anti-Rabbit IgG(H+L) Dylight 800	BS10026	Sheep anti-Mouse IgG(H+L) Dylight 549
BS10043	Goat Anti-Rat IgG(H+L) HRP	BS10038	Goat Anti-Rabbit IgG(H+L) Dylight 680	BS50251	Mouse Anti-Goat IgG (H+L) TRITC
BS10000	Goat Anti-Rat IgG(H+L) FITC	BS10034	Goat Anti-Rabbit IgG (H+L) Dylight 649	BS20509	Mouse Anti-Goat IgG (H+L) FITC
BS10009	Goat Anti-Rat IgG(H+L) Cy3	BS10029	Goat Anti-Rabbit IgG(H+L) Dylight 594	BS50551	Mouse Anti-Goat IgG (H+L) AP
BS10010	Goat Anti-Rat IgG(H+L) Biotin	BS10023	Goat Anti-Rabbit IgG(H+L) Dylight 549	BS30250	Rabbit Anti-Sheep IgG (H+L) TRITC
BS10007	Goat Anti-Rabbit IgG(H+L) Cy3	BS10017	Goat Anti-Rabbit IgG(H+L) Dylight 488	BS30950	Rabbit Anti-Sheep IgG (H+L) FITC
BS10005	Goat Anti-Rabbit IgG(H+L) AMCA	BS10012	Goat Anti-Rabbit IgG(H+L) Dylight 405	BS10030	Donkey Anti-Rabbit IgG(H+L) Dylight 594
BS10350	Goat Anti-Rabbit IgG (H+L)-HRP	BS10650	Goat Anti-Rabbit IgG (H+L) Biotin	BS10018	Donkey Anti-Rabbit IgG(H+L) Dylight 488
BS10003	Goat Anti-Rabbit IgG (H+L) Texas Red	BS10040	Goat Anti-Mouse IgG(H+L) Dylight 800	BS10028	Donkey Anti-Mouse IgG(H+L) Dylight 594
BS10250	Goat Anti-Rabbit IgG (H+L) TRITC	BS10037	Goat Anti-Mouse IgG(H+L) Dylight 680	BS10016	Donkey Anti-Mouse IgG(H+L) Dylight 488
BS10950	Goat Anti-Rabbit IgG (H+L) FITC	BS10027	Goat Anti-Mouse IgG(H+L) Dylight 594	BS10020	Donkey Anti-Goat IgG(H+L) Dylight 488
BS10006	Goat Anti-Mouse IgG(H+L) Cy3	BS10022	Goat Anti-Mouse IgG(H+L) Dylight 549	BS10903	Goat Anti-Human IgG (H+L) HRP
BS10004	Goat Anti-Mouse IgG(H+L) AMCA	BS10015	Goat Anti-Mouse IgG(H+L) Dylight 488	BS10995	Goat Anti-Human IgG (H+L) FITC
BS10002	Goat Anti-Mouse IgG (H+L) Texas Red	BS10011	Goat Anti-Mouse IgG(H+L) Dylight 405	BS10044	HRP Conjugated Stretavidin
BS11502	Goat Anti-Mouse IgG (H+L) TRITC	BS10033	Goat Anti-Mouse IgG (H+L) Dylight 649	Goat Anti-Rabbit	Goat Anti-Mouse

另有近萬種一抗可供選購，更多資訊請上：<https://www.bioworld.com/>





## Cytokine、重組蛋白特價中


Cat. No.	Product Name	Cat. No.	Product Name	Cat. No.	Product Name
BK0002	Recombinant Human Adiponectin	BK0066	Recombinant Human GM-CSF	BK0103	Recombinant Mouse IL-3
BK0003	Recombinant Mouse Adiponectin	BK0067	Recombinant Mouse GM-CSF	BK0106	Recombinant Human IL-5
BK0025	Recombinant Human EGF	BK0080	Recombinant Human IFN- $\gamma$	BK0107	Recombinant Mouse IL-5
BK0027	Recombinant Mouse EGF	BK0081	Recombinant Mouse IFN- $\gamma$	BK0108	Recombinant Rat IL-5
BK0028	Recombinant Rat EGF	BK0082	Recombinant Rat IFN- $\gamma$	BK0109	Recombinant Human IL-6
BK0053	Recombinant Human FGF2 (145aa)	BK0089	Recombinant Rat IL-10	BK0120	Recombinant Rat Leptin
BK0054	Recombinant Human FGF2 (146aa)	BK0091	Recombinant Human IL-15	BK0168	Recombinant Bovine TNF- $\alpha$
BK0055	Recombinant Human FGF2 (154aa)	BK0094	Recombinant Human IL-1 $\beta$	BK0170	Recombinant Mouse TNF- $\alpha$ , His
BK0056	Recombinant Bovine FGF2 (bFGF)	BK0095	Recombinant Mouse IL-1 $\beta$	BK0171	Recombinant Human TNF- $\alpha$
BK0057	Recombinant Mouse FGF2 (bFGF)	BK0097	Recombinant Rat IL-1 $\beta$	BK0184	Recombinant Human BDNF
BK0058	Recombinant Rat FGF2 (bFGF)	BK0098	Recombinant Human IL-21	BK0198	Recombinant Human EPO
BK0061	Recombinant Human G-CSF	BK0099	Recombinant Mouse IL-21	BK0208	Recombinant Human GDNF
BK0062	Recombinant Rat G-CSF	BK0102	Recombinant Human IL-3	BK0220	Recombinant Human HGF

找不到您要的品項？更多Cytokine、重組蛋白，請洽詢翰新國際！

或上 <https://www.bioworld.com/> 以獲得更多產品資訊

## 豐富的Plasmid/Vector產品線

Mammalian Expression Vector	Insect Expression Vector	Luciferase Reporter Vector
pCMV-Flag	pFastBac	pSTAT1-luc
pCMV-HA	pFastBac-His	pSTAT3-luc
pCMV-Myc		pNF $\kappa$ B-luc
pcDNA3.1-Flag	Cloning Vector	Lentiviral Expression Vector
pcDNA3.1-HA	pCMV	pLVX-Puro
pEGFP	pOTB7	Lentiviral Interference Vector
pEYFP		pLKO.1-Puro
pDsRed		
pmCherry		



感謝全世界用戶支持，PubMed引用已突破六千篇！

所有產品均通過嚴格品管測試，品質保證！產品優惠價格請洽當區業務或聯絡翰新：

