


只要A+B · 光速完成Western Blot配膠工序

不用再為了加錯APS, TEMED, Acrylamide而傷腦筋!!




不只快速凝膠，還附帶Gradient Gel特性!

普通12%膠

HyGel 1:1


改變A:B比例
可調整高、低分子量解析度

高分子量
不易分開



低分子量
容易跳海

VS



各分子量
均勻分布

A:B Ratio*	目標分子量
1:1	15~200KDa
2:3	>200KDa
3:2	<15KDa
7:3	<10KDa

*此比例僅供參考
比例不同凝膠速度也不同

< 訂購資訊 >

HyGel 1:1 Rapid Stacking Gel Forming Kit (上膠)	GB-SEBE250	2x250ML	\$3000
	GB-SEBE125	2x125ML	\$1800
HyGel 1:1 Rapid Resolving Gel Forming Kit (下膠)	GB-REBE500	2x500ML	\$9000
	GB-REBE250	2x250ML	\$5500

GOAL Bio Western Blot產品滿額好康大放送!!



凡購買WB系列產品滿15萬
即可抽iPhone 12等大獎!



產品編號	品名	包裝	售價
W-3400	HyBlock 1-min Blocking Buffer (PBST)	500 ml	\$2,800
W-3401	HyBlock 1-min Blocking Buffer (TBST)	500 ml	\$2,800
W-2105-8	Western Blot Stripping Buffer	500 ml	\$3,400
W-3408-0	HyECL Western Chemiluminescent Kit	2 x 250ml	\$5,800
W-3408-1	HyECL Western Chemiluminescent Kit	2 x 170ml	\$4,500
W-5XPD	5X SDS-PAGE Protein Sample Buffer	20 ml	\$1,800
W-2320-0	HyBradford 5X Protein Assay Reagent	500 ml	\$3,800
W-2320-7	7X Protein Gel HyStain Reagent	500 ml	\$3,800
W-7849	Hycell RIPA Buffer	500 ml	\$3,800
W-7850	Hytra Cell Protein Extraction Reagent	500 ml	\$6,800
W-7851	Hytra Tissue Protein Extraction Reagent	500 ml	\$7,800
W-7852	Hytra Mitochondrial Protein Extraction Kit	100 rxns	\$7,800
W-7883	Hytra Nuclear/Cytosol Protein Extraction Kit	100 rxns	\$7,800
HC100-007	100X Protease Inhibitor Cocktail	2ml	\$3,000
HC100-008	100X Phosphatase Inhibitor Cocktail	2ml	\$3,000
GB-RGB500	4x resolving gel buffer	500 ml	\$1,800
GB-SGB500	4x stacking gel buffer	500 ml	\$1,800
GB-SEBE250	HyGel 1:1 Rapid Stacking Gel Forming Kit	2 x 250 ml	\$3,000
GB-SEBE125	HyGel 1:1 Rapid Stacking Gel Forming Kit	2 x 125 ml	\$1,800
GB-REBE500	HyGel 1:1 Rapid Resolving Gel Forming Kit	2 x 500 ml	\$9,000
GB-REBE250	HyGel 1:1 Rapid Resolving Gel Forming Kit	2 x 250 ml	\$5,500
W100-032	10X TBST	5 x 1L	\$4,500
W100-003	10X PBST	5 x 1L	\$4,500
W100-014	10X TG-SDS (Running Buffer)	5 x 1L	\$4,500
W100-041	10X TG (Transfer Buffer)	5 x 1L	\$4,500
W100-020	20% SDS	500 ml	\$1,300
MB-63920-9	DimondTaq DNA Polymerase	500 rxns	\$21,000
AA0530-0100	Bovine Serum Albumin (BSA)	100g	\$7,800
GB-RN0200	GOALzol RNA Reagent	200 ml	\$5,900
PM2500	ExcelBand™ 3-color Regular Range Protein Marker	2 x 250 µl	\$3,200
MB-1166-8	HyFectin Transfection Reagent Kit	5 x 1 ml	\$45,000
MB-100BP	100bp DNA Ladder	5 x 500ul	\$5,000
MB-1KB	1kb DNA Ladder	5 x 500ul	\$5,000

