


無菌刻度 吸管



Non Pyrogenic
Non Cytotoxic
Non Haemolytic

- ◎ 五種體積可供選擇
- ◎ 單支無菌、紙背／膠膜包裝，方便易撕好取用
- ◎ 清楚且細微的刻度標示，一目了然
- ◎ 以不同顏色過濾棉，來區分不同體積
- ◎ Non pyrogenic、Non cytotoxic、Non haemolytic

規格	刻度	最大體積	包裝	特價
1 ml	1/100 ml	1.3 ml	500 's/pkg、4 pkg/cs	
2 ml	1/50 ml	2.4 ml	500 's/pkg、4 pkg/cs	
5 ml	1/10 ml	8 ml	250 's/pkg、4 pkg/cs	
10 ml	1/10 ml	13 ml	200 's/pkg、4 pkg/cs	
25 ml	1/5 ml	35 ml	150 's/pkg、4 pkg/cs	