



翰新國際

獨家代理羅氏

特惠期間至 7月31日

花東、輔大、桃園地區(不含長庚)

細胞分析實驗 細胞生理真實呈現

- Ready-to-use 試劑組方便直接使用，不需要額外花費時間進行調配。
- 簡易的操作流程，讓你輕鬆上手實驗。
- 設計快速，精準分析。

Product	Cellular Event							Detection Method	Key Benefit
	Apoptosis	Cell Proliferation	Necrosis	Flow Cytometry	Fluorescence Microscopy	Light Microscopy	Gel Electrophoresis		
Cell Death Detection ELISA ^{PLUS}	■	■					■	Easily discriminate between apoptosis and necrosis.	
Cell Proliferation Reagent WST-1		■					■	Achieve fivefold greater sensitivity compared to MTT or XTT.	
Cell Proliferation ELISAs BrdU		■					■	Accurately quantify cell proliferation without using radioactivity.	
ATP Bioluminescence Assay Kit HS II		■						Precisely measure even low concentrations of ATP.	
<i>In situ</i> Cell Death Detection Kits	■		■	■	■			Identify apoptosis at early stages at the molecular level using direct TUNEL detection.	
Annexin V-FLUOS Kit, Staining Kit, Annexin-V-Conjugates	■	■	■	■	■			Readily differentiate between apoptosis and necrosis by measuring membrane alterations.	
Apoptotic DNA-Ladder Kit	■	■				■		Easily distinguish apoptosis ("ladder") versus necrosis ("smear") using gel electrophoresis.	

更多資訊，歡迎參考 www.roche-applied-science.com/apoptosis

X-tremeGENE為羅氏最新一代超低毒性的轉染試劑，成功地轉染超過兩百種細胞株。

尋找各種細胞株適用的實驗操作步驟以及獲得目前使用者的使用心得及已發表的文獻，請上 www.x-tremegene.roche.com

30% off

X-tremeGene 給您超高速轉染效率及超低毒性 新一代轉染試劑，已成功轉染超過200株細胞株

Product Description	Cat. No.	Pack Size
Cell Cytotoxicity		
Cytotoxicity Det.Kit ^{PLUS} (LDH)	04 744 926 001	1 kit (400 tests in 96 wells)
Cell Proliferation		
Cell Proliferation ELISA, BrdU (colorimetric)	11 647 229 001	1 kit (1,000 tests)
Cell Proliferation Reagent WST-1	05 015 944 001	8 ml (800 tests)
	11 644 807 001	25 ml (2,500 tests)
5-Bromo-2'-deoxy-uridine Labeling and Detection Kit I	11 296 736 001	1 kit (100 tests)
Cell Death		
<i>In Situ</i> Cell Death Detection Kit, Fluorescein	11 684 795 910	1 kit (50 tests)
<i>In Situ</i> Cell Death Detection Kit, AP	11 684 809 910	1 kit (50 tests)
<i>In Situ</i> Cell Death Detection Kit, POD	11 684 817 910	1 kit (50 tests)
<i>In Situ</i> Cell Death Detection Kit, TMR red	12 156 792 910	1 kit (50 tests)
Cell Death Detection ELISA ^{PLUS}	11 774 425 001	1 kit (96 tests)
Transfection of mammalian cells		
X-tremeGENE HP DNA Transfection Reagent	06 366 244 001	0.4 ml
	06 366 236 001	1.0 ml
	06 366 546 001	5 x 1 ml
X-tremeGENE 9 DNA Transfection Reagent	06 365 779 001	0.4 ml
	06 365 787 001	1.0 ml
	06 365 809 001	5 x 1 ml

翰新國際有限公司

訂貨專線: (02) 2823-9957

e-mail: hycellbio@gmail.com



翰新國際

獨家代理羅氏

特惠期間至 7月31日

花東、輔大、桃園地區(不含長庚)

DIG 核酸標定系統 非放射線核標定與偵測

DIG 是羅氏專利. 更是專家

- **高靈敏度** - 可獲得與放射線相同或更高的靈敏度。
- **大幅節省操作時間** - 傳統放射線標定實驗流程需數日甚至數周，而DIG實驗流程僅需數分鐘至數小時。
- **保存更容易** - 無半衰期，僅需照一般核酸保存方式可達數年之久。
- **以建立標準實驗流程並有數千篇文獻** - DIG系統完全適用於膜雜交 (filter hybridization) 與原位雜交，我們並提供從探針合成至偵測的完整產品線

Product Description	Cat. No.	Pack Size
Labeling and Detection Kits for Northern Blots		
DIG Northern Starter Kit	12 039 672 910	1 kit (10 labeling reactions and detection of 10 blots)
Labeling and Detection Kits for Southern Blots		
DIG DNA Labeling and Detection Kit	11 093 657 910	1 kit (25 labeling reactions and detection of 50 blots)
DIG-High Prime DNA Labeling and Detection Starter Kit I	11 745 832 910	1 kit (12 labeling reactions and detection of 24 blots)
DIG-High Prime DNA Labeling and Detection Starter Kit II	11 585 614 910	1 kit (12 labeling reactions and detection of 24 blots)
Labeling and Detection Kits for Gel Shift Assays		
DIG Gel Shift Kit, 2 nd gen.	03 353 591 910	1 kit (20 labeling reactions and detection 20 blots)

Product Description	Cat. No.	Pack Size
Labeling Kits for RNA and DNA		
DIG RNA Labeling Kit (SP6/T7)	11 175 025 910	1 kit (2 x 10 labeling reactions)
PCR DIG Probe Synthesis Kit	11 636 090 910	1 kit (25 reactions of 50 µl final reaction volume)
Other Reagents to Compliment Northern, Southern and Gel Shift Assays		
Anti-Digoxigenin-AP, Fab fragments	11 093 274 910	150 U (200 µl)
DIG Easy Hyb	11 603 558 001	500 ml
	11 796 895 001	Granules for 6 x 100 ml
DIG Wash and Block Buffer Set	11 585 762 001	1 set (30 Blots)
CDP-Star, ready-to-use	12 041 677 001	2 x 50 ml
Lumi-Film Chemiluminescent Detection Film	11 666 916 001	100 films (18 x 24 cm)
	11 209 272 001	10 sheets (20 x 30 cm)
	11 209 299 001	20 Sheets (10 x 15 cm)
Nylon Membranes, positively charged	11 209 299 001	20 Sheets (10 x 15 cm)
	11 417 240 001	1 Roll (0.3 x 3 m)

翰新國際有限公司

訂貨專線: (02) 2823-9957

e-mail: hycellbio@gmail.com