



REAGEN LLC
 102 Commerce Dr. Unit 8
 Moorestown, NJ 08057 USA
 Tel: (856) 7270250; Fax: (856) 7270251
www.reagen.us

The Company

REAGEN LLC is an industry-leading biomedical enterprise focused on designing and producing a broad range of ELISA kits, including clinical diagnostic kits, food and feed safety kits, and veterinary disease kits. In support of our mission to enable cutting edge advances in the life sciences, we also manufacture microplate readers, plate washers, and related instruments and equipment. REAGEN is proud to provide customized programs to serve the specific needs of our clients, while maintaining our high standards and competitive prices.

Antibiotics and veterinary residues

Catalog	Product Name	Standards	Sensitivity	Quantity
RND99001	Chloramphenicol (CAP) ELISA Test Kit	0, 0.05, 0.15, 0.5, 1.5, 4.5ppb	0.05ppb	96T
RND99002	Streptomycin ELISA Test Kit	0, 0.5, 1, 2.5, 5, 10ppb	0.5ppb	96T
RND99003	Neomycin ELISA Test Kit	0, 1.5, 4.5, 13.5, 27, 81ppb	1.5ppb	96T
RND99004	Gentamicin (GM) ELISA Test Kit	0, 0.08, 0.32, 1.25, 5.0, 20ppb	0.08 ppb	96T
RND99005	Nitroimidazole ELISA Test Kit	0, 0.5, 1.5, 3, 6, 12ppb	0.5ppb	96T
RND99006	Tylosin ELISA Test Kit	0, 1, 2.5, 5, 10, 25ppb	1.0 ppb	96T
RND99007	Acrylamide-ES ELISA Test Kit	0, 2.5, 5, 10, 25, 50, 200ppb	2.5ppb	96T
RND99008	Malachite Green (MG) ELISA Test Kit	0,0.05,0.15,0.5,1.5,4.5 ppb	0.05ppb	96T
RND99009	Fluoroquinolone (FQS) ELISA Test Kit	0, 0.1, 0.25, 0.5, 1.0, 5.0ppb 0, 2.5, 7.5, 20, 50, 150ppb	0.1ppb(meat) 2.5ppb(honey)	96T
RND99010	Melamine ELISA Test Kit	0, 5, 20, 50, 100, 200ppb	5 ppb	96T
RND99011	Furazolidone (AOZ) ELISA Test Kit	0, 0.02, 0.06, 0.18, 0.54, 1.62ppb	0.02ppb	96T
RND99012	Furaltadone (AMAZ) ELISA Test Kit	0, 0.025, 0.1, 0.4, 1.6, 6.4ppb	0.05ppb	96T
RND99013	Clenbuterol ELISA Test Kit	0, 0.05, 0.15, 0.5, 1.5, 4.5ppb	0.05ppb	96T
RND99014	Ractopamine ELISA Test Kit	0, 0.5, 1.0, 5.0, 10, 25ppb	0.5 ppb	96T
RND99015	Sulfadiazine ELISA Test Kit	0, 1, 2, 5, 10, 25ppb	1ppb	96T
RND99016	Sulfamethazine ELISA Test Kit	0, 0.025, 0.5, 1, 5, 10ppb	0.025ppb	96T
RND99017	Sulfamethoxazole ELISA Test Kit	0, 0.5, 1.5, 5, 15, 50ppb	0.5ppb	96T
RND99018	Cigua-Check Fish Poison ELISA Kit	—	—	96T
RND99019	Sudan red ELISA Test Kit	0, 0.25, 0.5, 1.0, 3.5, 5ppb	0.25ppb	96T
RND99020	Ciprofloxacin ELISA Test Kit	0, 0.35, 1.0, 2.5, 7.5, 22.5ppb	0.35ppb	96T

RND99021	Ofloxacin ELISA Test Kit	0, 0.1, 0.25, 0.5, 1.0, 5ppb	0.1ppb	96T
RND99022	Chlorpromazine ELISA Test Kit	0, 0.04, 0.1, 0.25, 0.5, 1.5ppb	0.04ppb	96T
RND99023	Nitrofurazone (SEM) ELISA Test Kit	0, 0.025, 0.1, 0.4, 1.6, 6.4ppb	0.05ppb	96T
RND99024	Nitrofurantoin (AHD) ELISA Test Kit	0, 0.025, 0.1, 0.4, 1.6, 6.4ppb	0.025ppb	96T
RND99025	Tetracycline(TET) ELISA Test Kit	0, 0.05, 0.15, 0.3, 0.6, 1.5ppb	0.05ppb	96T
RND99026	Total Antibiotics in Milk Test Kit	—	—	96T
RND99027	Beta-lactam antibiotics ELISA Kit	0, 0.08, 0.2, 0.4, 0.8, 1.2ng/ml	0.04ppb	96T
RND99028	Reserpine ELISA Test Kit	0, 0.25, 0.5, 2.5, 5, 10ppb	0.25ppb	96T
RND99029	Florfenicol ELISA Test Kit	0, 0.15, 1, 5, 25, 100ppb	0.15ppb	96T
RND99030	Pasteurized verification ELISA Kit	—	—	96T
RND99031	The no-mel total milk protein ELISA Test Kit	—	—	96T
RND99032	Ivermectin ELISA Test Kit	0, 0.5, 2.0, 4.0, 8.0, 16.0 ppb	0.5ppb	96T
RND99033	Kanamycin ELISA Test Kit	0, 0.5, 1.5, 4.5, 13.5, 40.5ppb	0.5ppb	96T
RND99034	Ethoxyquin ELISA Test Kit	0, 0.1, 0.3, 0.9, 2.25, 4.5ppb	0.1ppb	96T
RND99035	Sulfonamides ELISA Test Kit	0, 0.5, 1.5, 5.0, 15, 50ng/ml	0.5ppb	96T
RND99036	Sulfaquinoxaline ELISA Test Kit	0, 1, 2.5, 5.0, 15, 50ng/ml	1.0ppb	96T
RND99037	Trimethoprim ELISA Test Kit	0,0.015,0.04,0.12,0.36,1.08 ppb	0.015ppb	96T
RND99038	Salbutamol ELISA Test Kit	0,0.05,0.15,0.5,1.5,4.5 ppb	0.05ppb	96T
RND99039	Cimaterol ELISA Test Kit	0, 0.25, 1.0, 5.0, 20, 50ppb	0.25ppb	96T
RND99040	Zilpaterol (ZIP) ELISA Test Kit	0, 0.015, 0.05, 0.1, 0.2, 0.5ppb	0.015ppb	96T
RND99041	Enrofloxacin ELISA Test Kit	0, 0.1, 0.25, 0.5, 1.0, 5.0ppb	0.1ppb	96T
RND99042	Flumequine ELISA Test Kit	0, 1, 2, 5, 15, 50 ng/ml	1.0ppb	96T
RND99043	Honey total antibiotic ELISA Kit	—	—	96T
RND99044	Crystal violet ELISA Test Kit	0,0.05,0.15,0.5,1.5,4.5 ppb	0.1ppb	96T
RND99045	Norfloxacin ELISA Test Kit	0, 1, 2, 5, 15, 50ng/ml	1ppb	96T
RND99046	Colistin ELISA Test Kit	0, 0.5, 1.5, 5, 15, 50ppb	0.5ppb	96T
RND99047	Aminoglycosides ELISA Test Kit	0, 0.5, 1, 2.5, 5, 10ppb	0.5ppb	96T
RND99048	Doxycycline ELISA Test Kit	0,0.05,0.15,0.5,1.5,4.5 ppb	0.05ppb	96T
RND99049	Oxytetracycline ELISA Test Kit	0,0.05,0.15,0.5,1.5,4.5 ppb	0.05ppb	96T
RND99050	Chlortetracycline ELISA Test Kit	0,0.05,0.15,0.5,1.5,4.5 ppb	0.05ppb	96T
RND99051	Olaquinox ELISA Test Kit	0, 1, 3, 9, 27, 81ppb	1ppb	96T
RND99052	Erythromycin ELISA Test Kit	0, 0.5, 1, 2, 4, 8 ppb	0.5ppb	96T
RND99053	Clonidine ELISA Test Kit	0,0.015,0.05,0.15,0.5,1.5,ppb	0.015ppb	96T
RND99054	Terbutaline ELISA Test Kit	0,0.05,0.15,0.5,1.5,4.5 ppb	0.05ppb	96T
RND99055	Cyproheptadine ELISA Test Kit	0,0.015,0.05,0.15,0.5,1.5,ppb	0.015ppb	96T
RND99056	Monensin ELISA Test Kit	0, 1, 2.5, 5.0, 7.5, 15.0 ppb	1ppb	96T
RND99057	Mabuterol ELISA Test Kit	0, 0.05, 0.15, 0.5, 1.5, 4.5ppb	0.05ppb	96T
RND99058	Lincomycin ELISA Test Kit	0, 0.2, 0.6, 1.8, 5.4, 16.2ppb	0.2ppb	96T
RND99059	Dexamethaone ELISA Test Kit	0, 0.1, 0.2, 0.4, 0.8, 1.6ppb	0.1ppb	96T

RND99060	Ceftiofur ELISA Test Kit	0, 2, 8, 16, 32, 64ppb	2ppb	96T
RND99061	Nicarbazin ELISA Test Kit	0, 1, 2, 4, 8, 16ppb	1ppb	96T
RND99062	Dimetridazole ELISA Test Kit	0, 0.5, 1.5, 3, 6, 12ppb	0.5ppb	96T
RND99063	Nalidixic acid ELISA Test Kit	0, 0.5, 1.5, 5, 15, 50ppb	0.5ppb	96T
RND99064	Tilmicosin ELISA Test Kit	0, 1, 2.5, 5, 10, 25ppb	1ppb	96T

The same kit have different Low Detection Limits for different samples. There are some examples as follows: *Sample Types And LOD For Antibiotics And Veterinary Residues Products(In ppb)*

Catalog	Product Name	Meat/Liver/Kidney	Milk	Honey	Feed
RND99001	Chloramphenicol (CAP) ELISA Test Kit	0.025	0.25	0.05	0.25
RND99004	Gentamicin (GM) ELISA Test Kit	2.5	2.5	6.25	6.25
RND99009	Fluoroquinolone (FQS) ELISA Test Kit	1.0	0.5	2.5	1.0
RND99011	Furazolidone (AOZ) ELISA Test Kit	0.04	0.04	0.04	0.1
RND99012	Furaltadone (AMAZ) ELISA Test Kit	0.05	0.05	0.05	0.1
RND99023	Nitrofurazone (SEM) ELISA Test Kit	0.05	0.05	0.05	0.1
RND99024	Nitrofurantoin (AHD) ELISA Test Kit	0.05	0.05	0.05	0.1

Mycotoxins

Catalog	Product Name	Standards	Sensitivity	Quantity
RNM98001	Aflatoxin M1 ELISA Test Kit	0, 0.05, 0.15, 0.5, 1.5, 5ppb 0, 0.05, 0.15, 0.3, 0.9, 2.7ppb	0.05ppb	96T
RNM98002	T2 Toxin ELISA Test Kit	0, 1, 5, 10, 25, 50 ppb	1.0ppb	96T
RNM98004	Deoxynivalenol ELISA Test Kit	0, 2.5, 5, 10, 25, 50 ppb	2.5ppb	96T
RNM98005	Zearalenone ELISA Test Kit	0, 0.1, 0.25, 0.5, 1.5, 4.5ng/ml	1.0ppb	96T
RNM98006	Total aflatoxin ELISA Test Kit	0, 0.05, 0.1, 0.2, 0.4, 0.8ng/ml	0.2ppb	96T
RNM98007	Aflatoxin B1 ELISA Test Kit	0, 0.05, 0.1, 0.2, 0.4, 0.8ng/ml	0.1ppb	96T
RNM98008	Ochratoxin A ELISA Test Kit	0, 0.015, 0.03, 0.06, 0.12, 0.24ppb	0.015ppb	96T
RNM98009	Fumonisin ELISA Test Kit	0, 1, 5, 10, 25, 50 ppb	1.0ppb	96T

Sample Types And LOD For Mycotoxins Products (In ppb)

Catalog	Product Name	Meat/Liver/Kidney	Milk	Feed
RNM98001	Aflatoxin M1 ELISA Test Kit	0.5	0.01	0.05
RNM98002	T2 Toxin ELISA Test Kit	3.5	10	3.5
RNM98004	Deoxynivalenol ELISA Test Kit	100	—	100
RNM98005	Zearalenone ELISA Test Kit	1.0	1.0	20
RNM98006	Total aflatoxin ELISA Test Kit	0.5	0.5	0.5
RNM98007	Aflatoxin B1 ELISA Test Kit	0.5/2.5	0.5	0.5/2.5
RNM98008	Ochratoxin A ELISA Test Kit	0.15	0.15	0.15
RNM98009	Fumonisin ELISA Test Kit	100	10	100

Algae Toxin

Catalog	Product Name	Standards	Sensitivity	Quantity
RNA97001	Cylindrospermopsin ELISA Kit	0,0.15, 0.4, 1.0, 2.0, 5.0ppb	0.04ppb	96T
RNA97002	Microcystins Tube Kit	0,0.15, 0.4, 1.0, 2.0, 5.0ppb	0.75ppb	40T
RNA97003	Microcystins ADDA ELISA Kit	0,0.15, 0.4, 1.0, 2.0, 5.0ppb	0.1ppb	96T
RNA97004	Microcystins DM ELISA Kit	0,0.15, 0.4, 1.0, 2.0, 5.0ppb	0.1ppb	96T
RNA97006	Brevetoxin (NSP) ELISA Test Kit	0, 0.010, 0.025, 0.05, 0.1, 0.25, 0.5, 2.0 ppb	0.01ppb	96T
RNA97007	Domoic Acid (ASP) ELISA Test Kit	0; 0.5; 1.0; 2.0; 5.0; 10.0 ppb	2.2ppb	96T
RNA97008	Okadaic Acid (DSP) ELISA Test Kit	0, 0.1, 0.2, 0.5, 1.0, 2.0, 5.0ppb	0.1ppb	96T
RNA97009	Saxitoxin (PSP) ELISA Test Kit	0, 0.02, 0.04, 0.08, 0.16, 0.32 ppb	0.02ppb	96T
RNA97010	Histamine Enzymatic Assay Kit	0, 1.0, 2.0, 4.0, 8.0, 12.0 ppm	1ppm	96T
RNA97011	Tetrodotoxin (TTX) ELISA Test Kit	0, 10, 20, 50, 200 ppb	10ppb	96T

Pesticide residues

Catalog	Product Name	Standards	Sensitivity	Quantity
RNP96001	2, 4-D ELISA Test Kit	0, 2.0, 10, 100 ppb	2.0ppb	96T
RNP96002	Acetochlor ELISA Test Kit	0, 0.1, 0.5, 2.5 ppb	42 ppt	96T
RNP96003	Alachlor ELISA Test Kit	0, 0.1, 0.25, 0.5, 1.0, 5.0ppb	1.0 ppb	96T
RNP96004	Atrazine ELISA Test Kit	0, 0.05, 0.1, 0.25, 1, 2.5, 5.0ppb	0.03 ppb	96T
RNP96006	Carbendazim ELISA Test Kit	0, 0.5, 1.5, 3.0, 6.0, 12ppb	0.5 ppb	96T
RNP96007	Cyclodienes ELISA Test Kit	0, 0.25, 0.5, 1.0, 2.5, 5, 10, 25 ppb	0.1ppb	96T
RNP96008	DDE/DDT ELISA Test Kit	0, 0.625, 1.25, 2.5, 5, 10, 25ppb	0.37ppb	96T
RNP96009	Diuron ELISA Test Kit	0, 0.03, 0.1, 0.3, 1.0, 3.0ppb	0.03 ppb	96T
RNP96010	Glyphosate ELISA Test Kit	0, 0.075, 0.2, 0.5, 1.0, 4.0ppb	0.05ppb	96T
RNP96011	Metolachlor ELISA Test Kit	0, 0.1, 0.25, 0.5, 1.0, 2.0, 5.0ppb	0.058ppb	96T
RNP96012	OP/C Screen Test Kit	Negative/Positive control	5ppb	20T / 96T
RNP96013	Pyrethroids ELISA Test Kit	0, 0.75, 2.5, 5.0, 15.0ppb	0.75 ppb	96T
RNP96014	Spinosyn ELISA Test Kit	0, 0.05, 0.125, 0.25, 0.5 ppb	0.05ppb	96T
RNP96015	Penoxsulam ELISA Test Kit	0, 0.25, 1.0, 10.0ppb	0.17ppb	96T
RNP96016	Fluridone ELISA Test Kit	0, 0.5, 2.0, 7.5, 15.0ppb	0.5ppb	96T
RNP96017	Trifluralin ELISA Test Kit	0, 0.25, 0.75, 1.5, 3.0, 7.5ppb	0.25ppb	96T

Industrial Pollutants

Catalog	Product Name	Standards	Sensitivity	Quantity
RNC95001	Bisphenol A ELISA Test Kit	0, 0.05, 0.2, 0.8, 2.5, 10 ppb	0.05ppb	96T
RNC95002	Coplanar PCBs ELISA Test Kit	0, 25, 50, 100, 250, 500, 1000ppt	25 ppt	96T
RNC95003	PCBs(Higher) Chlorinated ELISA Kit	0.25, 2.5, 25 ppb	0.25 ppb	100T
RNC95004	PBDE ELISA Test Kit	0, 25, 50, 100, 500, 1000 ppt	25 ppt	100T
RNC95005	Triclosan ELISA Test Kit	0.025, 0.1, 1.0 ppb 0,0.05,0.1,0.25,0.5,1.0,2.5ppb	0.02 ppb	100T/96T
RNC95006	Benzo[a]pyrene (BAP) ELISA TestKit	0, 0.25, 0.5, 1.0, 1.5, 2.5ppb	0.25ppb	96T

Hormone

Catalog	Product Name	Standards	Sensitivity	Quantity
RNH94001	17-Beta-Estradiol(E2) ELISA Kit	0, 3.0, 6.0, 9.0, 15.0, 25.0ppt	2.7ppt	96T/ 100T
RNH94002	Estrogens (E1+E2+E3) ELISA Ki	0, 0.05, 0.15, 0.5, 3.0 ppb	0.05ppb	96T
RNH94003	Estrone (E1) ELISA Test Kit	0, 0.05, 0.3, 0.8, 5.0ppb	0.05ppb	96T
RNH94004	Ethinylestradiol (EE2) ELISA Kit	0, 7.5, 15, 50, 100ppb	7.5ppb	96T
RNH94005	Estriol E3 ELISA Test Kit	0, 3, 6, 9, 15ppt	3ppt	96T
RNH94006	Methyltestosterone ELISA Test Kit	0, 0.1, 0.5, 2.5, 10, 50ppb	0.1ppb	96T
RNH94007	Testosterone ELISA Test Kit	0,0.05,0.15,0.45,1.35,4.05ppb	0.05ppb	96T
RNH94008	Boldenone ELISA Test Kit	0, 0.05, 0.15, 0.5, 1.5, 4.5ppb	0.05ppb	96T
RNH94009	Nandrolone ELISA Test Kit	0,0.015,0.05,0.15,0.5,1.5ppb	0.015ppb	96T
RNH94010	Trenbolone ELISA Test Kit	0, 0.05, 0.1, 0.25, 0.5, 1.0ppb	0.05ppb	96T
RNH94011	Diethylstilbestrol(DES) ELISA Kit	0, 0.15, 0.5, 1.5, 4.5, 9.0ppb	0.15	96T
RNH94013	Progesterone ELISA Test Kit	15.62, 31.25, 62.5, 125, 250,500ppt	15.62ppt	96T
RNH94014	Trenbolone ELISA Test Kit	0, 0.05, 0.1, 0.25, 0.5, 1.0ppb	0.05ppb	96T
RNH94015	Zeranol ELISA Test Kit	0, 0.1, 0.25, 0.5, 1.5, 4.5ppb	1ppb	96T

REAGEN

